

COST-EFFECTIVE EMCCD-BASED DETECTOR FOR TIME-RESOLVED BIOLOGICAL SAXS APPLICATIONS

V. Nagarkar, B. Singh, V. Gaysinskiy, L. Ovechkina, I. Shestakova, and S. Miller
Radiation Monitoring Devices, Inc.,
Watertown, MA 02472, USA

L. Guo, D. Gore, and T. Irving
BioCAT, Center for Synchrotron Radiation Research, Illinois Institute of Technology,
Chicago, IL 60616, USA

IWORID-8: July 2 – 6, 2006



Outline

- ❖ Background on time-resolved SAXS
- ❖ Detector requirements
- ❖ Detector concept and related developments
- ❖ Preliminary experiments at Argonne National Laboratory
- ❖ Results
- ❖ Discussion

Background on Time Resolved SAXS

- ❖ Time resolved X-ray diffraction of muscle
 - ★ Elementary force generation in muscle contraction
- ❖ Static and time-resolved scattering from macromolecules in solution
 - ★ Addresses the kinetics of “protein folding”
- ❖ Phase transitions in model membrane systems
 - ★ Understanding membrane fusion
 - ★ Biotechnological applications, especially in designing various drug delivery systems

Technical Challenges in X-Ray Data Acquisition in Partially Ordered and Disordered Systems

- ❖ Scattering from such systems is intrinsically weak
- ❖ Time scale of interest is ms or sub-ms level
- ❖ In solution, the weak scattering at higher angles is spread out isotropically – difficult to detect over the background
- ❖ In muscle, the fiber patterns from biological tissues are complex with the need to resolve closely spaced peaks over a wide intensity range
- ❖ Cover a wide solid-angle

Detector Requirements


❖ Scintillator

- ★ Bright
- ★ Fast decay time with no afterglow
- ★ High x-ray absorption efficiency
- ★ Excellent Spatial resolution
- ★ **Our Choice - New columnar CsI(Tl,Eu)**



❖ Readout Sensor

- ★ High pixel resolution
- ★ High frame rates
- ★ Low read and dark noise
- ★ Wide dynamic range
- ★ **Our Choice - New EMCCD Sensors**


The Detector Concept





RMD Columnar CsI(Tl)





Fast scintillator development and evaluation




RMD Columnar CsI(Tl, Eu)




RMD Columnar CsI(Tl, Eu)




RMD Columnar CsI(Tl,Eu)



Scintillator persistence: 30 ms/frame




Current CsI(Tl) Scintillator



Co-doped CsI(Tl,Eu)



Background on EMCCD



Advantages of EMCCD

- ❖ Low noise
 - Low read noise at high speeds through EM gain
- ❖ High speed readout
 - 30 fps to over 500 fps
- ❖ High sensitivity
 - Improved SNR for “light – starved” applications
- ❖ Adequate dynamic range

Custom Designed BI EMCCD




Custom Designed IXON 887 BI EMCCD



Detector Evaluation at the Argonne National Laboratory: Solution scattering of NtrC1C protein


- ❖ Used standard SAXS instrument on the BioCAT beamline 18-ID
- ❖ Camera length: 1.82 m
- ❖ X-ray wavelength: 0.103 nm
- ❖ Detectors used:
 - ★ Current EMCCD-based (512×512 pixels, 16×16μm² pixels)
 - ★ Existing state-of-the-art “Brandeis” detector (4k×7k pixels, 12×12μm² pixels)

Scattering Intensity (I) Versus Scattering Vector (Q) where $Q = (4\pi/\lambda) * \sin(\theta/2)$




Solution scattering of NtrC1C protein measured using the state-of-art detector and the new EMCCD detector with 3:1 Taper

Fiber Diffraction Pattern from the Embedded Muscle Specimen




**State-of-the-art
detector**

(30 ms data)



**EMCCD with 3:1
taper, no binning**

(33 ms data)



**EMCCD with 3:1
taper, 4×4 binning**

(11 ms data)

Z scale adjusted to account for higher EMCCD sensitivity

EMCCD Gain = 255


Summary

- ❖ Newly developed CsI(Tl,Eu)
 - ★ Preserves excellent emission properties of CsI(Tl)
 - ★ Enhances prompt emission by 20%
 - ★ Reduces afterglow by a factor of 30
- ❖ EMCCD are an excellent choice for SAXS application
 - ★ Fast readout rates with minimal read noise
 - ★ High frame rates of 500+ fps possible
- ❖ Large area fiberoptic coupled EMCCD detector development is currently underway at RMD
 - ★ Effective area $10.25 \times 10.25 \text{ cm}^2$
 - ★ High resolution of 1600×1600 pixels


Acknowledgements

- ❖ We thank the Public Health Service (NIH), DHHS for Grants R44-EB003382-01 and -02, and the Medical Sciences Div., DoE, for Grant DE-FG01-03ER83760, which provided partial support for parts of this work.


Quantum Accounting Diagram: Commercial VS EMCCD Based X-ray and γ -ray imaging systems



Quantum Accounting Diagram: Commercial VS EMCCD Based X-ray and γ -ray imaging systems



SNR at 10 MHz Readout Speed: Commercial VS EMCCD Based X-ray and γ -ray imaging systems




Imaging System Requirements

- ❖ Operate at the high global count rates ($\sim 10^7$ - 10^8 photons/s)
- ❖ High sensitivity (ability to detect single photons)
- ❖ High spatial resolution ($50\mu\text{m}$ - $100\mu\text{m}$)
- ❖ Large active area ($>10\text{cm}\times 10\text{cm}$)
- ❖ Read out with millisecond time resolution
- ❖ Wide dynamic range


EMCCD Specifications

<i>Parameter</i>	<i>Specification</i>
<i>CCD Chip</i>	E2V CCD887 BI
<i>Illumination</i>	Back Illuminated
<i>CCD Format</i>	512 x 512 Pixels
<i>Pixel Size</i>	16 μm Square
<i>Active Area (1:1 fiberoptic window)</i>	8.2 x 8.2 mm^2
<i>Active Area (3:1 fiberoptic window)</i>	24.6 x 24.6 mm^2
<i>Readout</i>	16 bits; 10MHz
<i>Frame rate (full resolution)</i>	32 fps
<i>Binning</i>	Flexible binning
<i>Frame rate binned mode</i>	Up to 520 fps
<i>On-chip gain</i>	1 to >1000x
<i>Operating temperature</i>	-30C
<i>Dark current</i>	1 e-/pixel/sec @-30C
<i>Read noise at 10 MHz</i>	30 e- (Unity gain)
<i>Read noise at 10 MHz</i>	<1 e- (Gain of 40)


High-speed test object




High Speed X-ray Imaging




27 fps




50 fps



80 fps



120 fps



224 fps